



NAUE GmbH

Sondermaschinen - und Filterbau

Am Rotböll 5
64331 Weiterstadt
www.naue.de

Tel. 0 61 50 / 97 38 - 0
Fax. 0 61 50 / 70 18
Info@naue.de

NAUE – Hose Connection (without any splits) for pharmacy and bio-technology

This hose connection represents the ideal connection to connect the hose with the connecting piece without any splits. Especially suitable for machines and installation sites.

It is one more possibility to optimize the quality assurance.

When repairing hoses of well-known manufacturers, we noted, that the so-called „hoses without any splits“ are obviously close, but often not “without any splits”.

During life of such a hose, plastics is flowing and splits are forming, which are for quality reasons not tolerable.

This mistake can be avoided by using the NAUE-Hose Connection.



Dimensions:

hose side: according to customers hose/specification

connecting side: sanitary fitting, clamp, welding end, etc.